Patient Information Leaflet


(Pancreatic drainage stent with Delivery System)

This leaflet has been written to help you understand stent. The doctor doing the procedure will discuss these in detail with you before they carry out the procedure.

* PRODUCT Information

Product Name	Pancreatic drainage stent with Delivery System
Brand Name	HANAROSTENT®
Manufacturer	M.I.Tech Co.,Ltd.
Model Name	HANAROSTENT® Plumber™
Intended use	This device is indicated for drainage of a pancreatic pseudocyst or a walled off necrosis through a transgastric or transduodenal approach.
Intended patient Population	Adult patients with pancreatic pseudocyst or walled off necrosis caused by pancreatic disease.
Expected lifetime	within 60 days ※ Since clinical data for patency after expected lifetime have not been verified, after the expected life time has elapsed, consult with your doctor to discuss follow up.

* What is PANCREATIC DRAINAGE STENT?

The Pancreatic drainage stent enables physicians to treat patients endoscopically instead of surgically. The Stent creates an anastomotic conduit between two lumens, enabling blockages and strictures to be bypassed and large fluid collections to be drained.

It have been primarily designed and developed to improve endoscopic treatment outcomes in the management of PFC, particularly walled-off necrosis (WON), and to facilitate better drainage of necrotic contents and minimize the risk of perforation and peritoneal or retroperitoneal leakage.

Patient Information Leaflet

(Pancreatic drainage stent with Delivery System)

* What are the risks involved?

Stent insertion is generally safe, but as with most medical treatment, there are some risks. These include: Some people get pain afterwards. This can be controlled with medicine if necessary. Occasionally a stent may slip out of position and the procedure will need repeating. Serious complications due to a stent insertion are rare. However, some side effects may require emergency treatment.

* **PREPARATION** for the procedure

Arrangements will be made for you to have a blood sample taken. This needs to be within two to three days before the procedure or if needed will be taken on the day of the test in the endoscopy unit. Please don't have anything to eat or drink for at least six hours before the procedure. If there are problems with your stomach emptying, this period may need to be longer but you will be told if this is necessary. You can take any normal medication swallowed with a small amount of water before the stent insertion. If you are diabetic or take warfarin, specific instructions regarding taking your medication will be given. If you are unsure what to do, please ask the procedure takes place at the endoscopy unit. On arrival an endoscopy nurse will meet you and carry out a health-check to make sure it is still safe to carry out the stent insertion on that day. The doctor will discuss with you before carrying out the procedure its intended benefits, risks of serious complications and any alternative treatment with you.

* How can you **CARE FOR YOURSELF**

- Rest when you feel tired. Getting enough sleep will help you recover.
- Avoid activities or exercises that use your belly muscles for 1 week or until your doctor says it is OK. For example, bicycle riding, jogging, weight lifting, or aerobic exercise.
- Don't lift or carry anything heavier than 4.5 kg (10 lb) for 3 days. As you feel ready, do a little more activity each day for the next 7 days after the procedure.
- Follow your doctor's directions for eating after the procedure.
- Drink plenty of fluids (unless your doctor tells you not to).
- If your doctor gave you a prescription medicine for pain, take it as prescribed.
- If you're not taking a prescription pain medicine, ask your doctor if you can take an over-the-counter medicine.
- Your doctor will tell you if and when you can restart your medicines. They will also give you instructions about taking any new medicines.
- If you take aspirin or some other blood thinner, be sure to talk to your doctor. They will tell you if and when to start taking this medicine again. Make sure that you understand exactly what your doctor wants you to do.
- **Follow-up care is a key part of your treatment and safety. Be sure to make and go to all appointments, and call your doctor or nurse call line if you are having problems.**
- **Please bring your patient implant card with you when taking MR imaging.**

* **SIDE EFFECTS**

The physician should be informed of any adverse events in patients who have had a stent implanted.

- Pain
- Bleeding
- Fever
- Rupture
- Infection

Patient Information Leaflet

(Pancreatic drainage stent with Delivery System)

- Tissue Necrosis
- Occlusion
- Stent misplacement or migration
- Perforation
- Edema
- Death
- Peritonitis
- Pancreatitis
- Jaundice
- Cholangitis
- Anesthesia complications
- Tissue ingrowth or overgrowth leading to difficult or a failure to remove stent
- Local infection at the implant site
- Sepsis (bacterial, endotoxin, or fungal)
- Persistent connection to fluid collection after removal (fistula)
- Cardia arrhythmia or arrest
- Partial or failed stent expansion, stent collapse
- Device failure, including failure to deliver the stent
- Abscess formation
- Aggravated walled off necrosis
- Adverse reaction to implant and/or delivery system (e.g., abdominal or back pain, nausea, infection, fever, chronic inflammation/foreign body reaction)
- Minor or excessive bleeding (requiring intervention)
- Leakage of fluid collection or bowel contents/peritonitis
- Tissue damage during stent implantation and/or removal
- Ulceration or erosion of mucosal or organ wall linings
- Pneumoperitoneum
- Stent fracture including damage of covering membrane

* CONTRAINDICATIONS

- All cardiovascular applications
- Hemodynamic instability
- Severe coagulopathy
- Cystic neoplasms.
- Pseudoaneurysms
- Duplication cysts
- Non-inflammatory fluid collections
- Patients with abnormal coagulation or who require ongoing complete anticoagulation at the time of implantation and post stent placement have an increased possibility of bleeding.
- Patients with altered anatomy that precludes the physician's ability to deliver the stent.
- Patients with intervening gastric varices or vessels within a one centimeter radius of the device needle.
- Allergy to metal (e.g. Nitinol, Nickel, Gold and Titanium)
- All others than indication for use

Patient Information Leaflet

(Pancreatic drainage stent with Delivery System)

* Incident notification

Notice that any serious incident that occurs in relation to the device should be reported to the manufacturer and to the regulatory authorities as below;

- Europe: European commission (<https://webgate.training.ec.europa.eu/eudamed-play/landing-page#/>)
- Australia: Therapeutic Goods Administration (www.tga.gov.au)
- Manufacturer: M.I.Tech Co., Ltd. (<http://www.mitech.co.kr>, mitech@mitech.co.kr)

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