

ST Hood

FUJIFILM

Distal End Cap

Small Caliber Tip

- Allows pocket creation method during ESD and POEM
- Tapered hood with water drains to support therapy

Maintains the endoscopic field of view during clinical procedures such as ESD

- Prevents the operative field from being disturbed by the mucous membrane
- Ensures a good operative field during endoscopic procedures such as EMR, ESD, POEM, etc



Small-caliber tip and transparent ST Hood for safer and more efficient ESD and POEM

Distal End Cap Features

Distal End Cap Features	ST Hood
• Providing a clear & stable view by keeping the distance between distal end & mucosal layer	✓
• Tapered design to facilitate tunneling into the submucosal layer	✓
• Side holes for better drainage of fluid and maintaining a clear endoscopic view	✓
• Sterile & single-use for better infection control	✓
• No guide ditch for easier insertion of accessories	✓



- **Wide inner diameter**
Good visibility can be ensured and a therapeutic device can be inserted smoothly.



- **Two drains**
The liquid inside the hood can be immediately poured out.

Ordering Information



Model	DH-28GR	DH-29CR	DH-30CR
Outer Diameter	11.5mm	12.4mm	13.8mm
Total Length	17.0mm	17.0mm	17.0mm
Distance from the Endoscope Tip	7.0mm	7.0mm	7.0mm
Inner Diameter of Distal End	8.0mm	8.0mm	8.0mm
Drain	2	2	2
Distal End Diameter of Compatible Endoscope	9.6 ~ 9.9mm (Gastroscope)	10.5 ~ 11.1mm (Slim Colonoscope)	11.7 ~ 12.8mm (Regular Colonoscope)
Compatible FUJIFILM Endoscope	EG-760Z EC-740T/M, L	EG-760CT	EC-760R-V/M, I, L EC-760ZP-V/M, L EC-720R/M, I, L

* Generic name: Single-use non-active equipment for natural orifice endoscopic treatment

* Package: 5pcs/pack