

# REN

## Balloon Dilatation Catheter

- Radial Balloon Dilatation
- Seamless Tip Design
- Visual Marker
- Radiopaque Marker





### Radial Balloon Dilatation

For preventing slippage from the papilla of Vater and providing withdrawal resistance



### Visual Marker

For easier positioning



### Seamless Tip Design

For smooth insertion into narrow parts



### Radiopaque Marker

For easier and precise balloon positioning under fluoroscopic guidance

## Structure



## Ordering Information

Product Code	Balloon Diameter	Balloon Length	Accepted GW	Sheath Length	Max Pressure	Working Channel
RN25-0320-18	3mm	20mm	0.025"	1,800mm	8atm	2.8mm
RN25-0430-18	4mm	30mm				
RN25-0630-18	6mm					
RN25-0830-18	8mm					
RN25-1240-18	10-11-12mm	40mm			7atm	3.2mm
RN25-1540-18	13-14-15mm					
RN25-1840-18	16-17-18mm					