


AXS_tome⁺

medw|ork[®]


SPHINCTEROTOME

- Micro-Tip⁺ design
- Self-locking mechanism (Handle brake)
- Multi-stage catheter



Micro-Tip+ design

- Smooth cannulation performance
- Minimizes tissue irritation and reduces risk of complications
- Optimized X-ray visibility
- Insulated cutting wire




Self-locking mechanism (Handle brake)



- Maintains the positioned flexion of the catheter tip when advancing the guidewire
- Easy access and precise tip control

Multi-stage catheter

- Maximum stability during intubation
- 3 Lumen



Ordering Information

Image	Product Code	Cutting Length	Distal Tip Length	Distal Tip Diameter	Sheath Length	Accepted GW	Working Channel
	PAP1-P8-25-35	25mm	5mm	4.2Fr	2,000mm	0.035"	2.8mm
	PAP1-P8-30-35	30mm	5mm	4.2Fr	2,000mm	0.035"	2.8mm