


HANAR  STENT®

# Biliary Flap (CCC)




“Malignant and Benign Biliary Obstruction”

- Fully covered stent with Anchoring Flaps

# HANAR STENT®

## Biliary Flap (CCC)

### Features



- 1 8 anchoring flaps for anti-migration
- 2 Gold radiopaque marker for precise stent placement

### Ordering Information

※ Sizing and availability varies by country

#### Endoscopic Application

Model	Stent (mm)			Delivery Device	
	Diameter	Usable Length	Total Length*	Length (mm)	Diameter (mm/fr)
BCG-06-040-180	6	20	40	1800	2.66/8
BCG-06-120-180		100	120	1800	2.66/8
BCG-08-040-180	8	20	40	1800	2.83/8.5
BCG-08-120-180		100	120	1800	2.83/8.5
BCG-10-040-180	10	20	40	1800	2.83/8.5
BCG-10-120-180		100	120	1800	2.83/8.5

\* Increment of 10mm

#### Percutaneous Application

Model	Stent(mm)			Delivery Device	
	Diameter	Usable Length	Total Length*	Length (mm)	Diameter (mm/fr)
BCG-06-040-060	6	20	40	600	2.66/8
BCG-06-120-060		100	120	600	2.66/8
BCG-08-040-060	8	20	40	600	2.83/8.5
BCG-08-120-060		100	120	600	2.83/8.5
BCG-10-040-060	10	20	40	600	2.83/8.5
BCG-10-120-060		100	120	600	2.83/8.5

\* Increment of 10mm