


CHOOSTENT[®] Esophagus (CCC)




“Esophageal Malignant Obstruction”


- Fully covered stent

CHOOSTENT[®] Esophagus (CCC)

Features



- 1 Dumbbell head on the both ends for anti-migration
- 2 Gold radiopaque marker for precise stent placement
- 3 Lasso on the both ends for stent repositioning



Segmented body structure to conform to peristaltic movement

Ordering Information

※ Sizing and availability varies by country

• NES

Model	Stent(mm)			Delivery Device	
	Diameter	Usable Length	Total Length*	Length (mm)	Diameter (mm/fr)
NES-18-060-070	24-18-24	20	60	700	6/18
NES-18-170-070		130	170	700	6/18
NES-20-060-070	26-20-26	20	60	700	6/18
NES-20-170-070		130	170	700	6/18
NES-22-060-070	28-22-28	20	60	700	8/24
NES-22-170-070		130	170	700	8/24
NES-24-060-070	30-24-30	20	60	700	8/24
NES-24-170-070		130	170	700	8/24

* Increment of 10mm