




EndoAll

Polypectomy Snare

- Braided wire
- 360° Rotation







Braided wire

Twisted multi-wire design, delivering enhanced grip and cutting performance for secure capture and safe resection.

360° Rotation


Smooth one-handed clockwise and counterclockwise rotation, allowing precise positioning and efficient lesion capture.



Ordering Information

A : Loop Diameter

* Wire Thickness



Product Image	Product Code	Type	Loop Diameter	Wire Thickness	Sheath O.D.	Working Length	Rotatable
 Oval	AG-5075-241023	Hybrid	10mm	0.24mm	2.4mm	2,300mm	Y
	AG-5075-241523	Hybrid	15mm	0.24mm	2.4mm	2,300mm	Y
	AG-5075-242523 **	Hot	25mm	0.24mm	2.4mm	2,300mm	Y
	AG-5075-243223	Hot	32mm	0.41mm	2.4mm	2,300mm	Y
 Shield	AG-5076-241023	Hybrid	10mm	0.24mm	2.4mm	2,300mm	Y
	AG-5076-241523	Hybrid	15mm	0.24mm	2.4mm	2,300mm	Y
 Round	AG-5078-241523	Hot	15mm	0.4-0.3-0.4mm*	2.4mm	2,300mm	Y
	AG-5078-242523	Hot	25mm	0.4-0.3-0.4mm*	2.4mm	2,300mm	Y

** For removal of EndoAll single-use hemoclip