

EndoAll

Biopsy Forceps

- Sharp cutting edges
- Wire hook design
- Sheath marking




Sharp cutting edges

Enabling clean and precise tissue sampling




Wire hook design

Providing stable cup control with reliable tissue capture



Sheath marking

Supporting accurate depth control and safe positioning



Ordering Information

Product Code	Type	Product Image	Cup Outer Diameter	Working Length	Min. Working Channel	Opening Width	Tissue Capacity	Cup Type	Swing	Rotation
AG-50121-2316	Standard Type		2.3mm	1,600mm	2.8mm	6.5mm	5.9mm ³	Oval	N	N
AG-50121-2323			2.3mm	2,300mm	2.8mm	6.5mm	5.9mm ³	Oval	N	N
AG-50121-1816			1.8mm	1,600mm	2.0mm	5.6mm	2.6mm ³	Oval	N	N
AG-5019-2316	Swing Type		2.3mm	1,600mm	2.8mm	6.5mm	5.9mm ³	Oval	Y	N
AG-5019-2323			2.3mm	2,300mm	2.8mm	6.5mm	5.9mm ³	Oval	Y	N
AG-5019-1816			1.8mm	1,600mm	2.0mm	5.6mm	2.6mm ³	Oval	Y	N
AG-5851-2316	Rotary Type		2.3mm	1,600mm	2.8mm	5.8mm	5.0mm ³	Oval	N	Y
AG-5851-2323			2.3mm	2,300mm	2.8mm	5.8mm	5.0mm ³	Oval	N	Y