


EndoAll



Injection Needle

- Sharp needle tip with dual-edge surface
- Ergonomic handle design
- Star-shaped inner catheter




Sharp needle tip with dual-edge surface

Ensuring high puncture force while preventing over-penetration




Ergonomic handle design

One-handed operation for easy and efficient procedures



Star-shaped inner catheter

Minimizing friction for smooth needle deployment and retraction



Ordering Information

Product Image	Product Code	Needle Gauge	Needle Length	Sheath O.D	Working Length
	AG-5311-2423-2304	23G (0.6mm)	4mm	2.4mm	2,300mm
	AG-5311-2423-2504	25G (0.5mm)	4mm	2.4mm	2,300mm